

Anwendungen Bio-Budget RNA-Isolierungskits

RNA-Isolierungssysteme							
Ausgangsmaterial	my Budget RNAmagic	my Budget RNA Mini Kit	my Budget Plant RNA Kit	my Budget micro RNA Kit	my Budget Blood RNA Mini Kit	my Budget DNA/RNA Mini Kit	my Budget Virus RNA Kit
Bakterien							
gram. pos.	✓	✓		✓		✓	
gram. pos.	✓	✓		✓		✓	
Blut							
Frisch/Vollblut	•	•			✓		
EDTA/Citrat-behandelt	•	•			✓		
Eukaryotische Zellen							
Monolayerzellen	✓	✓		✓		✓	
Suspensionszellen	✓	✓		✓		✓	
Zellkulturüberstand	✓	✓		✓		✓	
Formalin-/Paraffinmaterial							
Formalin	•	•					
Paraffin	•	•					
Gewebe/Biopsien							
Innere Organe	✓	✓		✓		✓	
Fettreiches Gewebe	✓	✓		✓		✓	
Hirnmateriale	✓	✓					
Würmer	✓	✓		✓		✓	
Fische	✓						
Pilze/Hefezellen							
Saccharomyces	✓		✓				
Pichia u.a.	✓		✓				
Pflanzenmaterial							
Gewebe	✓		✓				
Zellkultur	✓		✓				
Viren							
Serum	•				•		
Plasma	•				•		
Rachenabstriche							✓
Zellfreie Körperflüssigkeiten	•						
Dauer der Präparation	ca. 60 min.	ca. 15 - 40 min.	ca. 30 min.	ca. 15 - 40 min.	ca. 45 min.	ca. 15 - 40 min.	ca. 25 min.

✓ sehr gut geeignet

• bedingt geeignet